

Interactive FT

ACD/NMR Processor v12.01

Joe DiMartino
Advanced Chemistry Development, Inc.
Toronto, ON, Canada

Introduction

The instructions below describe how to pre-process your data. For more information on the manual processes that are not described below, please refer to the software manual for details.


Loading Raw Data

On the **File** menu, point to **Import**, and click **From 2D NMR Directory** or click  to open the **Import** dialog box.


Select a file format of **(Auto Detect)** and click drives and directories until you locate the **fid** file that you wish to process.

Interactive Fourier Transform

Click [Interactive FT](#) and the Fourier Transform dialog box will open (see screen shot below):



By selecting any of the headers (i.e. Linear Prediction, Window Function), this dialog box will expand out:



Fourier Transform

Horizontal Along t2

Apply FT along t2

Initial: 512 Final: 512

t2 FT Properties

Linear Prediction

Window Function

Vertical Along t1

Apply FT along t1

Initial: 128 Final: 512

t1 FT Properties (Dataset)

Linear Prediction

Backward

Points to Predict: 1 Base Points: 2 Coeff Count: 8

Forward

Points to Predict: 129 Base Points: 1 Coeff Count: 8

Window Function

Phasing Parameters

Auto Phasing

Mouse Phasing


Instant Preview

Preview


OK Cancel Help

If you would like to see an instant preview of the 2D NMR data, select **Instant Preview** and it will contain all the changes you made.

Note If you are not satisfied with how the contours appear, you can still do further processing with the [Interactive FT](#) dialog box.




In this case, select **Auto Phase** in the **Interactive FT** dialog box and it will automatically update the 2D NMR spectrum:



Click **OK** once you are pleased with the processing.

Note **Interactive FT** is also in the ACD/1D NMR Processor software and follows the same workflow as what was mentioned above.



Conclusion

The process described above is a basic overview of how to use Interactive FT for pre-processing 2D NMR data with ACD/2D NMR Processor. Many more processing options exist that were not described above. For a more in-depth view of any of the individual processing features, please consult the software manual. Please also visit www.acdlabs.com for more details.